

BonePlant is an affordable synthetic bone implant based on the latest developments in materials science.

Composition:


- 60 % HA (Hydroxyapatite)
- 40% TCP (Tricalcium Phosphate)
- The advantages of BonePlant:
 - excellent biocompatibility
 - perfect tissue integration
 - reliability




BonePlant Powder



Powder Bone, granule type, can be rolled down and can not be supplied by blood efficiently.



Teeth extraction. It makes big def-
fection and needs
bone grafting to
insert implants.



New bone is re-
generated within a
given period.