

EZ-Seal MTA Mineral Trioxide Aggregate



EZ-Seal is a calcium silicate cement-based root canal filling material which is used to induce restoration of pulp and periodontal tissue in root canal treatment at a dental clinic.

Advantages of EZ-Seal comparing to other root filling materials:

- aluminium was not added
- the product is heavy metal free
- EZ-Seal does not cause discoloration
- EZ-Seal's hardening time is short
- EZ-Seal is MTA which does not expand inside the root canal
- Inserting into the root canal is easy



ALUMINUM FREE
HEAVY METALS FREE
FAST SETTING TIME
HIGH BIOLOGICAL PERFORMANCE

Mineral Trioxide Aggregate
Regenerative Endodontic Repair Material

