

We offer a wide range of Medtronic (formerly Covidien, the leading medical device manufacturer in the world) suture materials for dental and oral surgery.

The needles are made from Surgalloy™ premium alloy (a new, high-tech nickel stainless steel), which provides high-quality performance and has two important properties: strength and flexibility.

NUCOAT™ technology - Surgical needles are coated with silicone. This reduces the initial penetration resistance force, thus ensuring greater precision and accuracy.

Reduces tissue pull on the needle body, minimizing tissue trauma for better wound healing. Improves the performance of multiple repeated tissue penetration, ensuring accuracy for the last stitch.

Wound closure products are available in absorbable, non-absorbable designs in different lengths and thicknesses, with a wide range of needle types.